



Miltenyi Biotec

Automated clinical-grade manufacturing of V γ 9V δ 2 T cells on the CliniMACS Prodigy[®]

Background

Gamma delta ($\gamma\delta$) T cells are a unique T lymphocyte subset crucial for immunosurveillance and tissue homeostasis.¹ The V γ and V δ chains of their T cell receptor (TCR) configuration enable MHC-independent recognition of diverse antigens.^{1,2} This feature sets them apart from $\alpha\beta$ T cells and makes them promising candidates for immunotherapy, particularly those using an allogeneic approach.¹⁻⁴

The overall presence of $\gamma\delta$ T cells in peripheral blood is low, with V γ 9V δ 2 T cells being the most common subset.² Using them in adoptive cell therapy requires establishing good manufacturing practice (GMP)-compliant methods to obtain sufficient cells for clinical doses. V γ 9V δ 2 T cells can be activated and expanded via exposure to natural or synthetic phosphoantigens, or aminobisphosphonate drugs such as zoledronic acid, in conjunction with cytokines.^{1,5}

Often, GMP-compliant protocols to expand V γ 9V δ 2 T cells involve multiple manual steps that risk contamination and yield insufficient cell numbers for clinical applications. Furthermore, while genetic engineering can enhance $\gamma\delta$ T cell cytotoxicity and targeting¹⁻³, their modification using common viral vector approaches remains challenging.

In this application note, we present data on the automated, clinical-grade manufacturing of V γ 9V δ 2 T cells using the CliniMACS Prodigy. The workflow developed consists of two processes (fig. 1). In the first process, TCR $\alpha\beta$ ⁺ and CD19⁺ cells are depleted from leukapheresis products to generate the starting material for subsequent expansion and optional transduction. In the second process, V γ 9V δ 2 T cells are preferentially activated and expanded from those depleted products using zoledronic acid. Additionally, activated V γ 9V δ 2 T cells can be successfully transduced within the same processing workflow using a baboon envelope pseudotyped (BaEV) lentiviral vector.




Figure 1: A two-process workflow for the generation of native or engineered $\gamma\delta$ T cells from allogeneic starting materials and their analysis, including in-process control and release testing.

Materials and methods

Depletion of $\alpha\beta$ T cells and B cells

Fresh, non-mobilized leukapheresis from healthy donors ($n = 3$) were depleted of $\alpha\beta$ T and B cells using the CliniMACS Prodigy LP-TCR α/β -19-45RA System. Depletion was performed in the single-use CliniMACS Prodigy Tubing Set (TS) 320 with CliniMACS[®] TCR α/β -Biotin, CliniMACS Anti-Biotin Reagent, and CliniMACS CD19 Reagent.

Expansion of $\gamma\delta$ T cells

The depleted cell fractions were transferred to the CliniMACS Prodigy TS 620 to proceed with Case 3 of CliniMACS Prodigy T Cell Transduction – Large Scale (TCT-LS) and a specific activity matrix for V γ 9V δ 2 T cell expansion. A total of 1×10^9 cells were seeded and cultured in stimulation medium consisting of TexMACS[™] GMP Medium containing 5% human AB serum, 140 U/mL MACS[®] GMP Recombinant Human IL-15, 500 U/mL MACS GMP Recombinant Human IL-2, and 2 μ g/mL zoledronic acid (pharmaceutical grade). As soon as the stimulation medium was consumed, medium without zoledronic acid was supplied, and the cells were maintained in culture for up to 14 days.

$\gamma\delta$ T cell transduction and VCN determination

Another set of depleted cell fractions ($n = 3$) were genetically modified via viral transduction. The cells were seeded and stimulated as described above using a different activity matrix that included automated transduction. On day 3 of stimulation, the cells were incubated for 24 h with BaEV lentiviral particles containing a specific chimeric antigen receptor (CAR) construct and 10 μ g/mL of MACS GMP Vectofusin[®]-1. Afterwards, the viral vector was washed out and the cells were expanded for up to 14 days. Vector copy number (VCN) was determined on day 14 using the MACS COPYcheck Kit, human, having isolated genomic DNA with the DNeasy[®] Blood & Tissue Kit (Qiagen).

QC testing by flow cytometry

The MACSQuant[®] Analyzer 10 was used to monitor T cell population growth and composition on days 0, 7, 10, and 14. Transduction efficiency was determined on days 7, 10, and 14 using the corresponding CAR Detection Reagent. The product table on the last page lists the antibodies and reagents used for analytical testing. For further details on the specific panels used, contact your local Technical Support.

Results

Efficient depletion of $\alpha\beta$ T and B cells to boost safety of allogeneic starting material for $\gamma\delta$ T cell expansion

$\alpha\beta$ T and B cells must be removed from allogeneic material for $\gamma\delta$ T cell therapy products to prevent graft-versus-host disease (GVHD) or B cell-triggered disorders.⁶ We performed this depletion before $\gamma\delta$ T cell expansion using an established automated process on the CliniMACS Prodigy, which resulted in a logarithmic depletion >4 for $\alpha\beta$ T cells and >3 for B cells (table 1). Consequently, $\gamma\delta$ T cells were partially enriched with a mean recovery of 87%. Given their low frequency in peripheral blood, however, that translates to 3.7% of the remaining cells. Other cells included monocytes and NK cells (table 1). On average, V γ 9V δ 2 T cells accounted for 72% of the $\gamma\delta$ T cell population.

Strong expansion and favorable phenotype of V γ 9V δ 2 T cells

Cell fractions depleted of $\alpha\beta$ T cells and B cells served as starting material for $\gamma\delta$ T cell expansion on the CliniMACS Prodigy. We used zoledronic acid added to the culture medium to stimulate V γ 9V δ 2 T cells. The combination with TexMACS GMP Medium, IL-2, IL-15, and optimized feeding conditions achieved expansion rates exceeding 200-fold. After 14 days in culture, 90% of the cells were $\gamma\delta$ T cells (98% of them were V γ 9V δ 2 T cells) with only low levels of $\alpha\beta$ T cells and B cells (table 1). Despite donor variation regarding expansion capacity, an average of 6.63×10^9 $\gamma\delta$ T cells were obtained by day 14, with up to 8.9×10^9 for one donor (fig. 2A). The expanded $\gamma\delta$ T cells were highly viable ($>97\%$) and exhibited a favorable central and effector memory phenotype (fig. 2B).

Parameter	Starting material depleted of $\alpha\beta$ T and B cells ($n = 3$)	Final cellular product after expansion ($n = 3$)
$\alpha\beta$ T cells (%)	0.007 \pm 0.004	0.130 \pm 0.180
$\gamma\delta$ T cells (%)	3.670 \pm 1.360	89.75 \pm 5.520
B cells (%)	0.005 \pm 0.003	0.007 \pm 0.005
Monocytes (%)	59.37 \pm 8.460	0.030 \pm 0.008
NK cells (%)	24.58 \pm 8.540	9.790 \pm 5.580
TCR $\alpha\beta^+$ cell log depletion	4.590 \pm 0.260	–
CD19 ⁺ cell log depletion	3.960 \pm 0.230	–

Table 1: Composition of starting material and final product of $\gamma\delta$ T cell manufacturing using the CliniMACS Prodigy Platform. Values represent mean and standard deviation.




Figure 2: Vγ9Vδ2 T cells are efficiently expanded on the CliniMACS Prodigy Instrument and exhibit a central and effector memory phenotype. By day 14 in culture, absolute cell numbers from the starting material of three healthy donors averaged 6.63×10^9 total viable $\gamma\delta$ T cells (A). The number of $\gamma\delta$ T cells with effector memory (CD27⁻CD45⁻) and central memory (CD27⁺CD45RA⁻) phenotype increased through expansion, reducing the proportion of naive (CD27⁺CD45RA⁺) and terminally differentiated (CD27⁻CD45RA⁺) cells.

BaEV lentiviral vector efficiently generates CAR $\gamma\delta$ T cells

In addition to establishing optimized expansion conditions for Vγ9Vδ2 T cells on the CliniMACS Prodigy, we evaluated the generation of CAR $\gamma\delta$ T cells by adapting the predefined culture conditions to include a viral transduction step. We used BaEV lentiviral particles to introduce a specific CAR construct. This choice was based on previous observations that BaEV lentiviral particles significantly increase transduction rates of $\gamma\delta$ T cells compared to the commonly used VSV-G pseudotype.⁷




Figure 3: CAR $\gamma\delta$ T cells are generated with high transduction efficiency and low VCN on the CliniMACS Prodigy.

The frequency of transduced $\gamma\delta$ T cells and NK cells expressing a specific CAR was evaluated on day 14 of culture (A). Likewise, the VCN of the final cell product was determined on day 14. Shown are values from three healthy donors.

By day 14 in culture, an average of 57% of the expanded $\gamma\delta$ T cells expressed the specific CAR (fig. 3A). Because many NK cells were still present when transduction took place on day 3, they were also efficiently transduced with BaEV lentiviral particles. Certain NK and CAR NK cells remained after the expansion period, contributing to the cytotoxic potential of the final cellular product (fig. 3A). After transduction, cells were expanded up to day 14. An average of 1.5×10^9 total $\gamma\delta$ T cells (among those 8.2×10^8 CAR $\gamma\delta$ T cells) were obtained with an average VCN of 1.4, suggesting a safe profile for the final cell product (fig. 3B).

Conclusions

This application note highlights the capabilities of the CliniMACS Prodigy in streamlining the manufacturing of V γ 9V δ 2 T cells to achieve clinically relevant doses for adoptive cell therapy, while enabling the effective depletion of $\alpha\beta$ T cells and B cells for safer use in allogeneic approaches. Moreover, the integration of viral transduction using BaEV lentiviral particles into the processing workflow enables the efficient generation of CAR $\gamma\delta$ T cells.

This automated and closed cell manufacturing workflow offers a promising avenue for the scalable and reproducible production of both unmodified and engineered V γ 9V δ 2 T cells for clinical applications, which represents a significant advancement for the development of novel cell-based immunotherapies.

References

- Handgretinger, R. and Shilbach, K. (2018) The potential role of $\gamma\delta$ T cells after allogeneic HCT for leukemia. *Blood* 131: 1063–1072.
- Hu, Y. *et al.* (2023) $\gamma\delta$ T cells: origin and fate, subsets, diseases and immunotherapy. *Signal Transduct. Target. Ther.* 8: 434.
- Saura-Esteller, J. *et al.* (2022) Gamma delta T-cell based cancer immunotherapy: past-present-future. *Front. Immunol.* 16: 915837.
- Yazdanifar, M. *et al.* (2020) $\gamma\delta$ T cells: the ideal tool for cancer immunotherapy. *Cells* 9: 1305.
- Bold, A. *et al.* (2023) An optimized cultivation method for future *in vivo* application of $\gamma\delta$ T cells. *Front. Immunol.* 14:1185564.
- Bethge, W. A. *et al.* (2021) Results of a multicenter phase I/II trial of TCR $\alpha\beta$ and CD19-depleted haploidentical hematopoietic stem cell transplantation for adult and pediatric patients. *Bone Marrow Transplant.* 57:423–430.
- Pinot, L. *et al.* (2025) Transduction of $\gamma\delta$ T cells with baboon envelope pseudotyped lentiviral vector encoding chimeric antigen receptors for translational and clinical applications. *Front. Immunol.* 16: 1548630.

Product table

Product	Order no.
CliniMACS Prodigy Instrument with LP-TCR $\alpha\beta$ -19-45RA Depletion application software and TCT-LS application software	200-075-301
CliniMACS Prodigy TS 320	Specific to regulatory variant required
CliniMACS TCR α/β -Biotin	
CliniMACS Anti-Biotin Reagent	
CliniMACS CD19 Reagent	
CliniMACS Prodigy TS 620	200-073-620
TexMACS GMP Medium, 2000 mL	170-076-306
MACS GMP Recombinant human IL-2	170-076-147
MACS GMP Recombinant human IL-15	170-076-114
MACS GMP Vectofusin-1	170-076-165
MACS COPYcheck Kit, human*	130-128-157
MACSQuant Analyzer 10	130-096-343
CD45 Antibody, anti-human, VioBlue®, REAfinity® (REA747)	130-110-637
TCR γ/δ Antibody, anti-human, APC, REAfinity (REA 591)	130-113-508
CD14 Antibody, anti-human, VioGreen™, REAfinity (REA599)	130-110-525
TCR α/β Antibody, anti-human, PE (BW242/412)	130-113-531
CD3 Antibody, anti-human, FITC, REAfinity (REA613)	130-113-138
CD56 Antibody, anti-human, PE-Vio® 615, REAfinity (REA196)	130-114-550
CD20 Antibody, anti-human, APC-Vio 770, REAfinity (REA780)	130-111-341
TCR V δ 2 Antibody, anti-human, VioBlue, REAfinity (REA771)	130-111-015
CD45RA Antibody, anti-human, VioGreen, REAfinity (REA562)	130-113-369
TCR V δ 1 Antibody, anti-human, FITC, REAfinity (REA173)	130-118-362
TCR V γ 9 Antibody, anti-human, APC-Vio 770, REAfinity (REA470)	130-128-218
CD27 Antibody, anti-human/mouse, PE, REAfinity (REA499)	130-113-640
7AAD Staining Solution	130-111-568
FcR Blocking Reagent	130-059-901
CAR Detection Reagents **	–

* Only compatible with the lentiviral vectors provided by Miltenyi Bioindustry and Miltenyi Biotec.

For more information on MACS antibodies, please visit [miltenyibiotec.com/antibodies](https://www.miltenyibiotec.com/antibodies)

** For information on CAR Detection Reagents, please visit [miltenyibiotec.com/CARdetection](https://www.miltenyibiotec.com/CARdetection)



Miltenyi Biotec

Miltenyi Biotec B.V. & Co. KG | Friedrich-Ebert-Straße 68 | 51429 Bergisch Gladbach | Germany | Phone +49 2204 8306-0 | Fax +49 2204 85197
macsde@miltenyi.com | www.miltenyibiotec.com

Miltenyi Biotec provides products and services worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Unless otherwise specifically indicated, Miltenyi Biotech products and services are for research use only and not for therapeutic or diagnostic use. For further regulatory notes, refer to [miltenyibiotec.com/regulatory-notes](https://www.miltenyibiotec.com/regulatory-notes). CliniMACS, CliniMACS Prodigy, MACS, MACSQuant, the Miltenyi Biotec logo, REAfinity, TexMACS, Vio, VioBlue, and VioGreen are registered trademarks or trademarks of Miltenyi Biotec B.V. & Co. KG and/or its affiliates in various countries worldwide. Vectofusin is a registered trademark of Genethon. All other trademarks mentioned in this document are the property of their respective owners and are used for identification purposes only. Copyright © 2025 Miltenyi Biotec and/or its affiliates. All rights reserved.